


## CSNK1D Antibody (R33607)

Catalog No.	Formulation	Size
R33607-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1453
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CSNK1D antibody is available for research use only.



Western blot testing of HEK293 cell lysate with CSNK1D antibody at 0.3ug/ml. Predicted molecular weight: ~47 kDa.

### Description

Additional name(s) for this target protein: Casein kinase 1, delta

### Application Notes

Optimal dilution of the CSNK1D antibody should be determined by the researcher.

## **Immunogen**

Amino acids DREERLRHSRNPATR were used as the immunogen for this CSNK1D antibody.

## **Storage**

Aliquot and store the CSNK1D antibody at -20oC.