


CSNK1E Antibody (R34962)

Catalog No.	Formulation	Size
R34962-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1454
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CSNK1E antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human CSNK1E tested with CSNK1E antibody [mock transfection in lane 1]. Predicted molecular weight: ~47kDa.

Description

Additional name(s) for this target protein: Casein kinase 1 epsilon

Application Notes

Optimal dilution of the CSNK1E antibody should be determined by the researcher.

Immunogen

Amino acids PASQTSVPFDHLGK were used as the immunogen for this CSNK1E antibody.

Storage

Aliquot and store the CSNK1E antibody at -20oC.