


## CSTB Antibody [Cystatin B] (R35488)

Catalog No.	Formulation	Size
R35488-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1476
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CSTB antibody is available for research use only.



Western blot testing of human lung lysate with CSTB antibody at 0.1ug/ml. Predicted molecular weight: ~11kDa, observed here at ~16kDa.

### Description

Additional name(s) for this target protein: CST6, cystatin B, stefin B

### Application Notes

Optimal dilution of the CSTB antibody should be determined by the researcher.

## Immunogen

Amino acids QTNKAKHDELTYP were used as the immunogen for this CSTB antibody.

## Storage

Aliquot and store the CSTB antibody at -20oC.