


CTnT Antibody (R35699)


Catalog No.	Formulation	Size
R35699-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Predicted Reactivity	Dog, Cat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7139
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CTnT antibody is available for research use only.



Western blot testing of rat heart lysate with CTnT antibody at 0.03ug/ml



Western blot of pig heart lysate with CTnT antibody at 0.01ug/ml.

Description

Additional name(s) for this target protein: Cardial troponin T, TNNT2, TnTC

Application Notes

Optimal dilution of the CTnT antibody should be determined by the researcher.

Immunogen

Amino acids RNRINDNQKVSKT were used as the immunogen for this CTnT antibody.

Storage

Aliquot and store the CTnT antibody at -20oC.