


CXCR5 Antibody / Burkitt lymphoma receptor 1 [clone CXCR5/13096] (V5994)


Catalog No.	Formulation	Size
V5994-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5994-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5994SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	CXCR5/13096
UniProt	P32302
Localization	Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml Western Blot : 2-4ug/ml
Limitations	This CXCR5/Burkitt lymphoma receptor 1 antibody is available for research use only.



SDS-PAGE Analysis of Purified CXCR5 Mouse Monoclonal Antibody (CXCR5/13096). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human adrenal gland stained with CXCR5 Mouse Monoclonal Antibody (CXCR5/13096). Inset: PBS instead of primary antibody; secondary only negative control.

Application Notes

Optimal dilution of the CXCR5/Burkitt lymphoma receptor 1 antibody should be determined by the researcher.

Immunogen

A recombinant fragment of human CXCR5 protein (exact sequence is proprietary) was used as the immunogen for the CXCR5/Burkitt lymphoma receptor 1 antibody.

Storage

CXCR5/Burkitt lymphoma receptor 1 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.