


CXCR5 Antibody / C-X-C motif chemokine receptor 5 [clone CXCR5/13100] (V5992)


Catalog No.	Formulation	Size
V5992-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5992-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5992SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)


Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	CXCR5/13100
UniProt	P32302
Localization	Cell membrane
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This CXCR5/C-X-C motif chemokine receptor 5 antibody is available for research use only.



SDS-PAGE Analysis of Purified CXCR5 Mouse Monoclonal Antibody (CXCR5/13100). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human brain stained with CXCR5 Mouse Monoclonal Antibody (CXCR5/13100). Inset: PBS instead of primary antibody; secondary only negative control.



Formalin-fixed, paraffin-embedded human adrenal gland stained with CXCR5 Mouse Monoclonal Antibody (CXCR5/13100). Inset: PBS instead of primary antibody; secondary only negative control.

Application Notes

Optimal dilution of the CXCR5/C-X-C motif chemokine receptor 5 antibody should be determined by the researcher.

Immunogen

A recombinant fragment of human CXCR5 protein (exact sequence is proprietary) was used as the immunogen for the CXCR5/C-X-C motif chemokine receptor 5 antibody.

Storage

CXCR5/C-X-C motif chemokine receptor 5 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.