


CXXC4 Antibody (R34144)

Catalog No.	Formulation	Size
R34144-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80319
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CXXC4 antibody is available for research use only.



IHC testing of FFPE human breast with CXXC4 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: CXXC finger 4

Application Notes

Optimal dilution of the CXXC4 antibody should be determined by the researcher.

Immunogen

Amino acids GTSLERTPVPSAEA were used as the immunogen for this CXXC4 antibody.

Storage

Aliquot and store the CXXC4 antibody at -20oC.