


CYB5R3 Antibody (R33271)


Catalog No.	Formulation	Size
R33271-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1727
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CYB5R3 antibody is available for research use only.



Western blot testing of human umbilical cord lysate with CYB5R3 antibody at 0.01ug/ml.
Predicted molecular weight (human): ~34/32/38kDa (isoforms 1/2/3).



IHC staining of FFPE human adrenal gland with CYB5R3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Cytochrome b5 reductase 3

Application Notes

Optimal dilution of the CYB5R3 antibody should be determined by the researcher.

Immunogen

Amino acids CLPNLDHVGHPTEP were used as the immunogen for this CYB5R3 antibody.

Storage

Aliquot and store the CYB5R3 antibody at -20°C.