


Cyb5r4 Antibody (R35548)

Catalog No.	Formulation	Size
R35548-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q68EJ0
Gene ID	266690
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Cyb5r4 antibody is available for research use only.



Western blot testing of rat spleen lysate with Cyb5r4 antibody at 0.1ug/ml. Predicted molecular weight (rat): ~59/53kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Cytochrome b5 reductase 4

Application Notes

Optimal dilution of the Cyb5r4 antibody should be determined by the researcher.

Immunogen

Amino acids RCQLEKLALREKRFD were used as the immunogen for this Cyb5r4 antibody.

Storage

Aliquot and store the Cyb5r4 antibody at -20oC.