


## Cyclin B2 Antibody (CCNB2) (F52193)

Catalog No.	Formulation	Size
F52193-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52193-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Primate
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	O95067
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This Cyclin B2 antibody is available for research use only.



### Description

Essential for the control of the cell cycle at the G2/M (mitosis) transition. [UniProt]

### Application Notes

Titration of the Cyclin B2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

A portion of amino acids 72-99 from the human protein was used as the immunogen for this Cyclin B2 antibody.

## **Storage**

Aliquot the Cyclin B2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.