


Cyclin B2 Antibody (CCNB2) (F52194)

Catalog No.	Formulation	Size
F52194-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52194-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O95067
Applications	Western Blot : 1:1000
Limitations	This Cyclin B2 antibody is available for research use only.



Description

Essential for the control of the cell cycle at the G2/M (mitosis) transition. [UniProt]

Application Notes

Titration of the Cyclin B2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 182-209 from the human protein was used as the immunogen for this Cyclin B2 antibody.

Storage

Aliquot the Cyclin B2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.