


## Cyclin B3 Antibody (CCNB3) (F52195)

Catalog No.	Formulation	Size
F52195-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52195-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q8WWL7
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This Cyclin B3 antibody is available for research use only.



### Description

CCNB3 is a mitotic cyclin that shares homology with cyclins of the A and B types. [OMIM]

### Application Notes

The stated application concentrations are suggested starting amounts. Titration of the Cyclin B3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 1311-1338 from the human protein was used as the immunogen for this Cyclin B3 antibody.

## Storage

Store at 4°C for up to one month. For long term, aliquot the Cyclin B3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.