


CYP2C8 Antibody (R34557)

Catalog No.	Formulation	Size
R34557-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1558
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CYP2C8 antibody is available for research use only.



Western blot testing of human liver lysate with CYP2C8 antibody at 0.1ug/ml. Predicted molecular weight: ~56kDa.

Description

Additional name(s) for this target protein: Cytochrome P450 2C8

Application Notes

Optimal dilution of the CYP2C8 antibody should be determined by the researcher.

Immunogen

Amino acids KIVALTRSYIREK were used as the immunogen for this CYP2C8 antibody.

Storage

Aliquot and store the CYP2C8 antibody at -20oC.