


Cystatin 3 Antibody (R34143)

Catalog No.	Formulation	Size
R34143-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	13010
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Cystatin 3 antibody is available for research use only.



Western blot testing of mouse testis lysate with Cystatin 3 antibody at 0.01ug/ml.
Predicted molecular weight: ~16 kDa.

Description

Additional name(s) for this target protein: Cystatin C; Cst3

Application Notes

Optimal dilution of the Cystatin 3 antibody should be determined by the researcher.

Immunogen

Amino acids GTHSLTKFSCKNA were used as the immunogen for this Cystatin 3 antibody.

Storage

Aliquot and store the Cystatin 3 antibody at -20oC.