


Cystatin C Antibody (R36096)

Catalog No.	Formulation	Size
R36096-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1471
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cystatin C antibody is available for research use only.



Western blot testing of human testis lysate with Cystatin C antibody at 0.1ug/ml.
Predicted molecular weight: ~16kDa.

Description

Additional name(s) for this target protein: CST3, ARMD11; CST3

Application Notes

Optimal dilution of the Cystatin C antibody should be determined by the researcher.

Immunogen

Amino acids SVEEEGVRRALD were used as the immunogen for this Cystatin C antibody.

Storage

Aliquot and store the Cystatin C antibody at -20oC.