

Cystatin F Antibody (R35176)

Catalog No.	Formulation	Size
R35176-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8530
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Cystatin F antibody is available for research use only.



Western blot testing of NSO lysate with Cystatin F antibody at 2ug/ml. Predicted molecular weight: ~17kDa.

Description

Additional name(s) for this target protein: Cystatin 7, CST7, CMAP

Application Notes

Optimal dilution of the Cystatin F antibody should be determined by the researcher.

Immunogen

Amino acids PWLQHFEVPVLRCH were used as the immunogen for this Cystatin F antibody.

Storage

Aliquot and store the Cystatin F antibody at -20oC.