


Cyt19 Antibody (R33106)

Catalog No.	Formulation	Size
R33106-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57412
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cyt19 antibody is available for research use only.



IHC testing of FFPE human kidney with Cyt19 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: AS3MT; arsenic (+3 oxidation state) methyltransferase

Application Notes

Optimal dilution of the Cyt19 antibody should be determined by the researcher.

Immunogen

Amino acids AALRDAEIQKDVQ were used as the immunogen for this Cyt19 antibody.

Storage

Aliquot and store the Cyt19 antibody at -20oC.