


CYTIP Antibody (R35969)

Catalog No.	Formulation	Size
R35969-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9595
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CYTIP antibody is available for research use only.



Western blot testing of MOLT4 with CYTIP antibody at 1ug/ml. Predicted molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Cytohesin-1 interacting protein

Application Notes

Optimal dilution of the CYTIP antibody should be determined by the researcher.

Immunogen

Amino acids SLQRLQHSNGN were used as the immunogen for this CYTIP antibody.

Storage

Aliquot and store the CYTIP antibody at -20oC.