


## Cytochrome b5 reductase 3 Antibody / CYB5R3 [clone 28C75] (RQ8486)

Catalog No.	Formulation	Size
RQ8486	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

<b>Availability</b>	1-3 days
<b>Species Reactivity</b>	Human, Rat
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	28C75
<b>Purity</b>	Affinity chromatography
<b>UniProt</b>	P00387
<b>Applications</b>	Western Blot : 1:500-1:2000
<b>Limitations</b>	This Cytochrome b5 reductase 3 antibody is available for research use only.



Western blot testing of human A431 cell lysate with Cytochrome b5 reductase 3 antibody. Predicted molecular weight ~34 kDa.

### Description

Cytochrome b5 reductase 3 catalyzes the reduction of two molecules of cytochrome b5 using NADH as the electron donor. [UniProt]

### Application Notes

Optimal dilution of the Cytochrome b5 reductase 3 antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide from human CYB5R3 protein was used as the immunogen for the Cytochrome b5 reductase 3 antibody.

## **Storage**

Store the Cytochrome b5 reductase 3 antibody at -20oC.