


## Cytochrome b5 reductase 3 Antibody / CYB5R3 [clone 28C75] (RQ8486)

Catalog No.	Formulation	Size
RQ8486	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human, Rat
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28C75
Purity	Affinity chromatography
UniProt	P00387
Applications	Western Blot : 1:500-1:2000
Limitations	This Cytochrome b5 reductase 3 antibody is available for research use only.



Western blot testing of human A431 cell lysate with Cytochrome b5 reductase 3 antibody. Predicted molecular weight ~34 kDa.

### Description

Cytochrome b5 reductase 3 catalyzes the reduction of two molecules of cytochrome b5 using NADH as the electron donor. [UniProt]

### Application Notes

Optimal dilution of the Cytochrome b5 reductase 3 antibody should be determined by the researcher.

### Immunogen

A synthetic peptide from human CYB5R3 protein was used as the immunogen for the Cytochrome b5 reductase 3 antibody.

## **Storage**

Store the Cytochrome b5 reductase 3 antibody at -20oC.