


Cytokeratin 19 Antibody Biotin Conjugate (R34300BTN)

Catalog No.	Formulation	Size
R34300BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3880
Applications	Western Blot : 0.3-1ug/ml Enzyme ImmunoAssay : suitable ELISA (peptide) LOD : 1:128000
Limitations	This Cytokeratin 19 antibody is available for research use only.



Western blot testing of human colon cancer lysate with biotinylated Cytokeratin 19 antibody at 0.3ug/ml. Predicted/observed molecular weight ~43/40-43kDa.

Description

Additional name(s) for this target protein: Keratin 19, CK19; KRT19

Application Notes

Optimal dilution of the Cytokeratin 19 antibody should be determined by the researcher.

Immunogen

Amino acids HTEQLQMSRSEVTD were used as the immunogen for this Cytokeratin 19 antibody.

Storage

Aliquot and store the Cytokeratin 19 antibody at -20oC.