


## Cytokeratin 19 Antibody / KRT19 [clone RM364] (R20383)

Catalog No.	Formulation	Size
R20383-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	RM364
<b>Purity</b>	Protein A purified from animal origin-free supernatant
<b>UniProt</b>	P08727
<b>Localization</b>	Cytoplasmic, cell membrane, secreted
<b>Applications</b>	Immunohistochemistry (FFPE) : 1:500-1:1000 Western Blot : 1:1000-1:2000
<b>Limitations</b>	This recombinant Cytokeratin 19 antibody is available for research use only.



IHC staining of FFPE human breast cancer tissue with recombinant Cytokeratin 19 antibody at 1:1000.



Western blot testing of human MCF7 cell lysate with recombinant Cytokeratin 19 antibody. Predicted molecular weight ~43 kDa.

## Description

This antibody reacts human CK-19 (Cytokeratin-19).

## Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Cytokeratin 19 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A peptide corresponding to the C-terminus of human CK-19 (Cytokeratin-19) was used as the immunogen for the recombinant Cytokeratin 19 antibody.

## Storage

Store the recombinant Cytokeratin 19 antibody at -20oC.