


Cytokeratin 19 Antibody (R34300)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34300-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 3880 |
| Applications | Western Blot : 0.03-0.1ug/ml Enzyme ImmunoAssay : used as capture/detect for recombinant protein ELISA (peptide) LOD : 1:128000 |
| Limitations | This Cytokeratin 19 antibody is available for research use only. |



Western blot testing of human colon cancer lysate with Cytokeratin 19 antibody at 0.03ug/ml. Predicted/observed molecular weight ~43/40-43kDa.

Description

Additional name(s) for this target protein: Keratin 19, CK19; KRT19

Application Notes

Optimal dilution of the Cytokeratin 19 antibody should be determined by the researcher.

Immunogen

Amino acids HTEQLQMSRSEVTD were used as the immunogen for this Cytokeratin 19 antibody.

Storage

Aliquot and store the Cytokeratin 19 antibody at -20oC.