


Cytokeratin 19 Antibody (R34300)

Catalog No.	Formulation	Size
R34300-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3880
Applications	Western Blot : 0.03-0.1ug/ml Enzyme ImmunoAssay : used as capture/detect for recombinant protein ELISA (peptide) LOD : 1:128000
Limitations	This Cytokeratin 19 antibody is available for research use only.



Western blot testing of human colon cancer lysate with Cytokeratin 19 antibody at 0.03ug/ml. Predicted/observed molecular weight ~43/40-43kDa.

Description

Additional name(s) for this target protein: Keratin 19, CK19; KRT19

Application Notes

Optimal dilution of the Cytokeratin 19 antibody should be determined by the researcher.

Immunogen

Amino acids HTEQLQMSRSEVTD were used as the immunogen for this Cytokeratin 19 antibody.

Storage

Aliquot and store the Cytokeratin 19 antibody at -20oC.