


D4-GDI Antibody (R36157)

Catalog No.	Formulation	Size
R36157-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	397
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This D4-GDI antibody is available for research use only.



Western blot testing of Jurkat lysate with D4-GDI antibody at 0.1ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: D4 GDP dissociation inhibitor, Rho GDP-dissociation inhibitor 2, ARHGDIB

Application Notes

Optimal dilution of the D4-GDI antibody should be determined by the researcher.

Immunogen

Amino acids TEKAPEPHVE were used as the immunogen for this D4-GDI antibody.

Storage

Aliquot and store the D4-GDI antibody at -20oC.