


D4-GDI Antibody (R36157)

Catalog No.	Formulation	Size
R36157-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	397
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This D4-GDI antibody is available for research use only.



Western blot testing of Jurkat lysate with D4-GDI antibody at 0.1ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: D4 GDP dissociation inhibitor, Rho GDP-dissociation inhibitor 2, ARHGDIB

Application Notes

Optimal dilution of the D4-GDI antibody should be determined by the researcher.

Immunogen

Amino acids TEKAPEPHVE were used as the immunogen for this D4-GDI antibody.

Storage

Aliquot and store the D4-GDI antibody at -20oC.