


DACH2 Antibody (R36272)

Catalog No.	Formulation	Size
R36272-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	117154
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This DACH2 antibody is available for research use only.



Western blot testing of mouse NIH3T3 cell lysate with DACH2 antibody at 0.03ug/ml.
Predicted/observed molecular weight (mouse): ~68kDa.

Description

Additional name(s) for this target protein: Dachshund 2

Application Notes

Optimal dilution of the DACH2 antibody should be determined by the researcher.

Immunogen

Amino acids TRKQAVNSSRPGR were used as the immunogen for this DACH2 antibody.

Storage

Aliquot and store the DACH2 antibody at -20oC.