


DAO Antibody (R35718)

Catalog No.	Formulation	Size
R35718-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1610
Applications	Western Blot: 0.03-0.01ug/ml ELISA (peptide) LOD: 1:16000
Limitations	This DAO antibody is available for research use only.



Description

Additional name(s) for this target protein: D-amino-acid oxidase

Application Notes

Optimal dilution of the DAO antibody should be determined by the researcher.

Immunogen

Amino acids RPQIRLEREQLRT were used as the immunogen for this DAO antibody. **Storage** Aliquot and store the DAO antibody at -20oC.