


DAO Antibody (R35718)

Catalog No.	Formulation	Size
R35718-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1610
Applications	Western Blot : 0.03-0.01ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DAO antibody is available for research use only.



Western blot testing of human cerebellum lysate with DAO antibody at 0.03ug/ml.
Predicted molecular weight: ~39kDa.

Description

Additional name(s) for this target protein: D-amino-acid oxidase

Application Notes

Optimal dilution of the DAO antibody should be determined by the researcher.

Immunogen

Amino acids RPQIRLEREQLRT were used as the immunogen for this DAO antibody.

Storage

Aliquot and store the DAO antibody at -20oC.