


## DARC Antibody Biotin Conjugate (R34265BTN)

Catalog No.	Formulation	Size
R34265BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2532
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : suitable Flow Cytometry : suitable ELISA (peptide) LOD : 1:16000
Limitations	This DARC antibody is available for research use only.



Western blot testing of human liver lysate with biotinylated DARC antibody at 0.1ug/ml.  
Predicted molecular weight: ~36 kDa.

### Description

Additional name(s) for this target protein: Duffy blood group / FY

### Application Notes

Optimal dilution of the DARC antibody should be determined by the researcher.

## **Immunogen**

Amino acids HRAELSPSTENSSQLDFED were used as the immunogen for this DARC antibody.

## **Storage**

Aliquot and store the DARC antibody at -20oC.