


DARC Antibody Biotin Conjugate (R34265BTN)

Catalog No.	Formulation	Size
R34265BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2532
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : suitable Flow Cytometry : suitable ELISA (peptide) LOD : 1:16000
Limitations	This DARC antibody is available for research use only.



Western blot testing of human liver lysate with biotinylated DARC antibody at 0.1ug/ml.
Predicted molecular weight: ~36 kDa.

Description

Additional name(s) for this target protein: Duffy blood group / FY

Application Notes

Optimal dilution of the DARC antibody should be determined by the researcher.

Immunogen

Amino acids HRAELSPSTENSSQLDFED were used as the immunogen for this DARC antibody.

Storage

Aliquot and store the DARC antibody at -20oC.