


# Dbf4 Antibody (R33506)

Catalog No.	Formulation	Size
R33506-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27214
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Dbf4 antibody is available for research use only.



# **Description**

Additional name(s) for this target protein: Dumbbell former 4

## **Application Notes**

Optimal dilution of the Dbf4 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DKKRTEYLPAHEDRT were used as the immunogen for this Dbf4 antibody.

Storage Aliquot and store the Dbf4 antibody at -20oC.						