


Dbf4 Antibody (R33506)

Catalog No.	Formulation	Size
R33506-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27214
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Dbf4 antibody is available for research use only.



Western blot testing of mouse testis lysate with Dbf4 antibody at 1ug/ml. Predicted molecular weight: ~74kDa.

Description

Additional name(s) for this target protein: Dumbbell former 4

Application Notes

Optimal dilution of the Dbf4 antibody should be determined by the researcher.

Immunogen

Amino acids DKKRTEYLP AHEDRT were used as the immunogen for this Dbf4 antibody.

Storage

Aliquot and store the Dbf4 antibody at -20oC.