


DDX6 Antibody (R33128)

Catalog No.	Formulation	Size
R33128-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1656
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DDX6 antibody is available for research use only.



Western blot testing of Daudi [A], Jurkat [B] and HepG2 [C] lysates with DDX6 antibody at 0.3ug/ml. Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: DEAD-box protien 6

Application Notes

Optimal dilution of the DDX6 antibody should be determined by the researcher.

Immunogen

Amino acids AEYHSEPVEDEKP were used as the immunogen for this DDX6 antibody.

Storage

Aliquot and store the DDX6 antibody at -20oC.