


DDX6 Antibody (R33128)

Catalog No.	Formulation	Size
R33128-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1656
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DDX6 antibody is available for research use only.



Western blot testing of Daudi [A], Jurkat [B] and HepG2 [C] lysates with DDX6 antibody at 0.3ug/ml. Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: DEAD-box protien 6

Application Notes

Optimal dilution of the DDX6 antibody should be determined by the researcher.

Immunogen

Amino acids AEYHSEPVEDEKP were used as the immunogen for this DDX6 antibody. **Storage** Aliquot and store the DDX6 antibody at -20oC.