


## Decorin Antibody (R34485)

Catalog No.	Formulation	Size
R34485-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1634
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Decorin antibody is available for research use only.



Western blot testing of human kidney lysate with Decorin antibody at 0.5ug/ml. Predicted molecular weight ~40 kDa (non-glycanated).

## Description

Additional name(s) for this target protein: DCN, CSCD, DSPG2, PG40, PGII, PGS2, SLRR1B; DCN

## Application Notes

Optimal dilution of the Decorin antibody should be determined by the researcher.

## Immunogen

Amino acids KISRVDAAASKGLNN were used as the immunogen for this Decorin antibody.

## **Storage**

Aliquot and store the Decorin antibody at -20oC.