


Decorin Antibody (R34485)

Catalog No.	Formulation	Size
R34485-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1634
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Decorin antibody is available for research use only.



Decorin Antibody Human Kidney WB. Western blot testing of human kidney lysate with Decorin antibody at 0.5ug/ml. Predicted molecular weight ~40 kDa (non-glycanated).

Description

Additional name(s) for this target protein: DCN, CSCD, DSPG2, PG40, PGII, PGS2, SLRR1B; DCN

This antibody can be compared with our [Decorin Antibody](#) (clone DCN/3521) for consistent detection of DCN across extracellular matrix and proteoglycan biology studies.

Application Notes

Optimal dilution of the Decorin antibody should be determined by the researcher.

Immunogen

Amino acids KISRVDAAASKGLNN were used as the immunogen for this Decorin antibody.

Storage

Aliquot and store the Decorin antibody at -20oC.