

DEF6 Antibody (R33147)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33147-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Dog, Mouse, Pig, Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 50619 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000 |
| Limitations | This DEF6 antibody is available for research use only. |



Western blot testing of human thymus lysate with DEF6 antibody at 1ug/ml. Predicted molecular weight: ~74kDa.

Description

Additional name(s) for this target protein: Differentially expressed in FDCP 6

Application Notes

Optimal dilution of the DEF6 antibody should be determined by the researcher.

Immunogen

Amino acids ALRKELLKSIWYA were used as the immunogen for this DEF6 antibody.

Storage

Aliquot and store the DEF6 antibody at -20oC.