


## Defensin 1/3 Antibody (R33429)


Catalog No.	Formulation	Size
R33429-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1667; 1668
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Defensin 1/3 antibody is available for research use only.



IHC staining of FFPE human thymus with Defensin 1/3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of human spleen lysate with Defensin 1/3 antibody at 2ug/ml. The expected ~10kDa band and the additional ~20kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: DEFA1; DEFA3

## Application Notes

Optimal dilution of the Defensin 1/3 antibody should be determined by the researcher.

## Immunogen

Amino acids DESLAPKHPGSRKN were used as the immunogen for this Defensin 1/3 antibody.

## Storage

Aliquot and store the Defensin 1/3 antibody at -20oC.