


Delta-Sarcoglycan Antibody Biotin Conjugate (R35676BTN)

Catalog No.	Formulation	Size
R35676BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow, Dog, Mouse, Pig
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6444
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:16000
Limitations	This delta-Sarcoglycan antibody is available for research use only.



Western blot testing of human heart lysate with biotinylated delta-Sarcoglycan antibody at 0.3ug/ml. Predicted molecular weight: ~35kDa.

Description

Additional name(s) for this target protein: SGCD; sarcoglycan, delta (35kDa dystrophin-associated glycoprotein); 35DAG

Application Notes

Optimal dilution of the delta-Sarcoglycan antibody should be determined by the researcher.

1. This delta-Sarcoglycan antibody will detect isoforms 1 and 3.

Immunogen

Amino acids RLPHGSYTPTGTR were used as the immunogen for this delta-Sarcoglycan antibody. This sequence is specific to isoforms 1 and 3.

Storage

Aliquot and store the delta-Sarcoglycan antibody at -20oC.