


Dennd4c Antibody (R35610)

Catalog No.	Formulation	Size
R35610-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	329877
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Dennd4c antibody is available for research use only.



Western blot testing of mouse spleen lysate with Dennd4c antibody at 0.5ug/ml.
Predicted molecular weight: ~211kDa.

Description

Additional name(s) for this target protein: DENN/MADD domain containing 4C

Application Notes

Optimal dilution of the Dennd4c antibody should be determined by the researcher.

Immunogen

Amino acids RFKPQSPYRAYKDRS were used as the immunogen for this Dennd4c antibody.

Storage

Aliquot and store the Dennd4c antibody at -20oC.