


DENTT Antibody (R35078)

Catalog No.	Formulation	Size
R35078-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64061
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This DENTT antibody is available for research use only.



IHC testing of FFPE human human testis with DENTT antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Differentially Expressed Nucleolar TGF beta1 Target, TSPY-like 2; TSPYL2

Application Notes

Optimal dilution of the DENTT antibody should be determined by the researcher.

Immunogen

Amino acids PNGWANPGKRGKTG were used as the immunogen for this DENTT antibody.

Storage

Aliquot and store the DENTT antibody at -20oC.