


## DHX9 Antibody (R36156)

Catalog No.	Formulation	Size
R36156-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1660
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DHX9 antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with DHX9 antibody at 1ug/ml. Predicted molecular weight: ~140kDa.

## Description

Additional name(s) for this target protein: DEAH (Asp-Glu-Ala-His) box polypeptide 9

## Application Notes

Optimal dilution of the DHX9 antibody should be determined by the researcher.

## Immunogen

Amino acids TEGRNALHKSSVNC were used as the immunogen for this DHX9 antibody.

## Storage

Aliquot and store the DHX9 antibody at -20oC.