


DHX9 Antibody (R36156)

Catalog No.	Formulation	Size
R36156-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1660
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This DHX9 antibody is available for research use only.



Description

Additional name(s) for this target protein: DEAH (Asp-Glu-Ala-His) box polypeptide 9

Application Notes

Optimal dilution of the DHX9 antibody should be determined by the researcher.

Immunogen

Amino acids TEGRNALIHKSSVNC were used as the immunogen for this DHX9 antibody. **Storage** Aliquot and store the DHX9 antibody at -20oC.