


DISP1 Antibody (R35881)

Catalog No.	Formulation	Size
R35881-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84976
Applications	Immunofluorescence : suitable ELISA (peptide) LOD : 1:32000
Limitations	This DISP1 antibody is available for research use only.



Immunofluorescence testing of U87 cells with DISP1 antibody at 2.5ug/ml

Description

Additional name(s) for this target protein: Dispatched homolog 1, DISPA

Application Notes

Optimal dilution of the DISP1 antibody should be determined by the researcher.

Immunogen

Amino acids SELSGESLLIKTL were used as the immunogen for this DISP1 antibody.

Storage

Aliquot and store the DISP1 antibody at -20oC.