


# DKC1 Antibody (R34622)

Catalog No.	Formulation	Size
R34622-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1736
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DKC1 antibody is available for research use only.



IHC testing of FFPE human spleen with DKC1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### **Description**

Additional name(s) for this target protein: Dyskeratosis congenita 1

## **Application Notes**

Optimal dilution of the DKC1 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids KRKRESESESDETPP were used as the immunogen for this DKC1 antibody. **Storage** Aliquot and store the DKC1 antibody at -20oC.