


DKC1 Antibody (R34622)

Catalog No.	Formulation	Size
R34622-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1736
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DKC1 antibody is available for research use only.



IHC testing of FFPE human spleen with DKC1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Dyskeratosis congenita 1

Application Notes

Optimal dilution of the DKC1 antibody should be determined by the researcher.

Immunogen

Amino acids KRKRESESESEDETPP were used as the immunogen for this DKC1 antibody.

Storage

Aliquot and store the DKC1 antibody at -20oC.