


DKK1 Antibody / Dickkopf-related protein 1 [clone FHG-4] (RQ4799)

Catalog No.	Formulation	Size
RQ4799	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	FHG-4
Purity	Affinity purified
UniProt	O94907
Applications	Western Blot : 1:1000-1:2000
Limitations	This DKK1 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with DKK1 antibody. Expected molecular weight: 26-40 kDa depending on glycosylation level.

Description

DKKs play an important role in vertebrate development, where they locally inhibit Wnt regulated processes such as antero-posterior axial patterning, limb development, somitogenesis and eye formation. [UniProt]

Application Notes

Optimal dilution of the DKK1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human DKK1 was used as the immunogen for the DKK1 antibody.

Storage

Store the DKK1 antibody at -20oC.