


DKK1 Antibody (R33211)


Catalog No.	Formulation	Size
R33211-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22943
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This DKK1 antibody is available for research use only.



Western blot testing of human bone marrow lysate with DKK1 antibody at 1ug/ml.
Predicted molecular weight: 26-40 kDa depending on glycosylation level.



IHC testing of FFPE human placenta tissue with DKK1 antibody at 10ug/ml. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

Description

Additional name(s) for this target protein: Dickkopf related protein-1

Application Notes

Optimal dilution of the DKK1 antibody should be determined by the researcher.

Immunogen

Amino acids CDHHQASNSSLHT were used as the immunogen for this DKK1 antibody.

Storage

Aliquot and store the DKK1 antibody at -20oC.