

## DKK4 Antibody (R35730)

Catalog No.	Formulation	Size
R35730-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27121
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This DKK4 antibody is available for research use only.



Western blot testing of Jurkat lysate with DKK4 antibody at 1ug/ml. Detection of the 30-32kDa form of the protein is observed. (1)

## Description

Additional name(s) for this target protein: Dickkopf homolog 4

## Application Notes

Optimal dilution of the DKK4 antibody should be determined by the researcher.

1. DKK4 is secreted as 40, 30-32, and 15-17 kDa proteins (OMIM: 605417).

## **Immunogen**

Amino acids RQHARLRVCQKIEKL were used as the immunogen for this DKK4 antibody.

## **Storage**

Aliquot and store the DKK4 antibody at -20oC.